

			"KEIKO"	"YAE"	"Cherry"
	Manufacturer		Kyoto Kagaku	Kyoto Kagaku	Kyoto Kagaku
	MPN		11405-000	11404-000	11401-000
	Web Link		<a href="https://us.vwr.com/store/product/22423983/kyoto-kagaku-patient-care-simulator-keiko">https://us.vwr.com/store/product/22423983/kyoto-kagaku-patient-care-simulator-keiko</a>	<a href="https://us.vwr.com/store/product/22423956/kyoto-kagaku-patient-care-simulator-yae">https://us.vwr.com/store/product/22423956/kyoto-kagaku-patient-care-simulator-yae</a>	<a href="https://us.vwr.com/store/product/22423714/kyoto-kagaku-patient-care-simulator-cherry">https://us.vwr.com/store/product/22423714/kyoto-kagaku-patient-care-simulator-cherry</a>
	VWR Part Number		470229-844	470229-842	470229-816
	Fidelity Rating		Low	Low	Low
	Function		Patient Care	Patient Care	Patient Care
Appearance	Size		149 cm approx.	152 cm approx.	62 in
	Weight		15 kg approx.	17 kg approx.	28.2 lb.
	Exchangeable Male/Female Genitalia		Female (Male Genitalia Sold Separately)	Female (Male Genitalia Sold Separately)	Male, Female
	Wig For Hair Care Exercises And Surgical Draping		Basic Hair Care	Basic Hair Care	X
	Replicable Mask For Age Appearance Change		X	Elder Mask & Wig (Sold Separately)	X
	Soft, Realistic Hands, Feet, Fingers, And Toes		X	X	X
	Carrying Case		BOX	BOX	BOX
Patient Care Tasks	Patient Hygiene	Bed Baths	X	X	X
		Hair Care Washing And Combing	X	X	X
		Clothing Change	X	X	X
		Denture Placement And Removal	X	X	X
		Movable Jaw With Tongue	X	X	X
		Oral And Nasal Hygiene	X	X	X
		Eye Irrigation	X	X	X
		Pressure Sore Care Stage 2 To 4, And 1 Unable To Judge	(heel, back and sacrum bone only)	(heel, back and sacrum bone only)	(heel, back and sacrum bone only)
		Suppository Insertion	X	X	X
		Anal Hemorrhoids	X	X	X
	Patient Positioning	Patient Transfer Techniques Capable	X	X	X
		Bends At Waist	X	X	X
		Sits Unassisted	X	X	X
		Joins With Realistic Range Of Motion	X	X	X
		Stoma Care	X	X	X
	Women's Health	Sculpted Stomas	Normal	Normal	Normal
		Pap Smears	X	X	X
		Realistic Vagina And Cervix	X	X	X
		Vaginal Douche	X	X	X
	Nutrition	Episiotomy 2Nd Degree Midline	X	X	X
		Nasogastric Tube Placement (For Tube Feeding)	(Tube Insertion Only)	(Tube Insertion Only)	(Tube Insertion Only)
Mouth, Nose And Tracheotomy Openings For Nasogastric Lavage And Gavage		X	X	X	
Gastrostomy Care		X	Not Available	X	
PEG Care		X	Not Available	X	
Bandaging Activity	TPN Care	X	X	X	
Physical Assessment	Eye and Pupillary Response	Manually Set Eyelids	X	X	X
		Blinking Eyelids	X	X	X
		Realistic Eyes In Eye Sockets For Ophthalmic Exercises	X	X	X
		Light Sensitive Pupils	X	X	X
Puncture / Injection	Injection Sites		X	X	X
		Injection Locations	Arm, thigh, buttock	Arm, thigh, buttock	Arm, Subcutaneous
	IV Insertion		X	X	Arm
		Iv Location	X	X	X
Catheter Openings		X	X	X	
Catheterization / Others	Urinary Catheterization		X	X	X
		Detachable And Removable Internal Tanks	X	X	X
		New Modular Tricuspid Valve Permits Male Or Female Catheterization With Soft Silicone Catheters	X	X	X
		Realistic Urethral Passage And Bladder For Catheterization Exercises	X	X	X
		Elevating Pillow Increases Bladder Pressure For Increased Flow During Catheterization (Male Or Female)	X	X	X
	Enema Administration	X	X	X	
Stomach Pumping		Not Available	X	X	
Airway Management and Intubation	Airway Intubation	Oral And Nasal Intubation	X	X	X
		Separate Left And Right Lungs	Not Available	X	Not Available
		Esophageal Intubation	X	X	X
		Suctioning Capabilities	X	X	X
		Endotracheal Tubes (Insertion, Securing)	Not Available	X	Not Available
		Spontaneous Breathing W/O External Compressor	X	X	X
		Supraglottic Airways	X	X	X
	Use Of A Real Ventilation Mask	X	X	X	
	Airway Suction	X	X	X	
	Tracheostomy Care	(Tube Insertion Only)	(Tube Insertion Only)	(Tube Insertion Only)	
Trach Opening For Placement	X	X	X		
Needle Chest Decompression	X	X	X		
CPR and Defib	CPR	Fully Articulated Head, And Neck	X	X	X
		Palpable And Visual Landmarks	X	X	X
		Jaw Thrust	X	X	X
		Mouth To Mouth CPR	X	X	X
		Mouth To Nose CPR	X	X	X
		Anatomical Landmarks	X	X	X
		Chest Compression Training	X	X	Not Available
Optional Trauma Materials	Wound Care Assessment Module	X	X	X	