



SONOtrain™ Cleaning


Ultrasound scan gel residue can be removed with a damp cloth or under running water.

After drying the model, sprinkle on the talcum powder provided. Rub it over the entire surface, the block is now ready for storage.



SONOtrain™ Construction

- Each individual unit is completely self-contained and can be used out of the box
- Units are stackable to save space
- Replacement blocks and all accessories are separately available



The vein model enables liquids to flow realistically through all the vessels shown.



9000768 09/2016

SONOtrain™

The ultrasound training set for medical simulation



THE ULTRASOUND TRAINING SET FOR MEDICAL SIMULATION

The material of SONOtrain™ simulates the feel of real soft tissue with realistic texture and echogenicity on ultrasound images.

The SONOtrain™ series is an excellent teaching tool for basic ultrasound training and for improving dexterity and hand-eye coordination skills.

It has been developed in Germany, is easy to use and a very cost effective teaching tool. Each training block can be used multiple times thanks to a „self-sealing“ material.

SONOtrain™ series for clinical ultrasound skill training:

- Realistic texture and echogenicity
- Train punctures and injections
- Perform tumor biopsies and cyst aspiration
- Diagnose gallbladder pathologies
- Detect foreign bodies in soft tissue

Ultrasound devices and needles not included in delivery content.




Item no. P120 (1019637)

SONOtrain™ Vein Model
 Ultrasound block with three blood vessels with diameters of 4, 8 and 15 mm with adjustable fluid flow. Punctures and injections possible, and self-closing injection channels.





Item no. P121 (1019636)

SONOtrain™ Foreign Body Model
 Ultrasound block with 6 different foreign bodies for insertion: Pellets, shard of glass, nail, projectile, loose chippings and splinters of wood can be clearly detected and recognized.




Item no. P122 (1019638)

SONOtrain™ Gallbladder Model
 Includes three gallbladders, each with different pathologies: Gallstones (diameter: 8 and 10 mm), thickened gallbladder wall and biliary sludge deposit. The gallstones and the biliary sludge deposit can move around.

Item no. P124 (1019634)

SONOtrain™ Breast Model with Cysts
 Realistic reproduction of a breast made from ultrasound material with two cysts for realistic punctures and aspirations.




Item no. P125 (1019635)

SONOtrain™ Breast Model with Tumours
 Realistic reproduction of a breast made from ultrasound material with three tumours for realistic punctures and biopsies.