



3B SCIENTIFIC® BIRTH STAGES TRAINER

Carry out realistic prenatal cervical examinations

Determine the labour progress using 6 inserts with various levels of cervical dilation and effacing.



The following levels and birth stages are shown:

- Cervix closed (0.5 cm), level -2
- Cervix 2 cm dilated, level -2
- Cervix 3 cm dilated, half effaced, level -2
- Cervix 5 cm dilated, effaced, level -1

- Cervix 7 cm dilated, effaced, level -1
- Cervix 9 cm dilated, effaced, level 0



Effective Training

The outside of the birthing stages trainers are neutral and do not indicate the birthing stage simulated inside. Only a number on the back indicates the stage.



Set contents:

- 6 inserts
- Safety plate with 4 screws (for optional insertion into the birthing simulator P90)
- Lubricant
- Carry bag



Each insert measures $13 \times 11.5 \times 13$ cm. Particularly durable and easy-to-clean material was used especially for the requirements of medical training.

TIP:

For more thorough and realistic training, the birth stage inserts can be used in combination with the birthing stimulator P90 BASIC and PRO.



Order list:

- > P94 / 1020628 3B Scientific® birth stages trainer, set
- > XP90-015 / 1020608 Lubricant (2 x 250 ml) – can also be ordered individually

May also be combined with:

- ▶ P90B / 1020332 3B Scientific® P90 Birthing Simulator BASIC
- > P90P / 1020333 3B Scientific® P90 Birthing Simulator PRO

10/2016 9000769

For further products relevant to extensive obstetric care training visit us online at **3bscientific.com**.

