



EPIDURAL AND SPINAL INJECTION TRAINER

- Extremely realistic haptic feedback
- High-quality material, needles don't leave a trace
- Practice loss-of-resistance and hanging-drop technique on standard or geriatric patient
- Fluid-filled spinal canal with realistic outflow rate



Main functions:

- Epidural anaesthesia using the loss-of-resistance and the hanging-drop technique
- Spinal anaesthesia with realistic resistance of the dura and arachnoid mater with or without a cannula
- Fluid-filled spinal canal with realistic outflow rate thanks to the ability to adjust the excess pressure
- An epidural catheter can be inserted into the epidural space
- Closed water system
- Easy to clean



> NEW

Geriatric Option

Simulates the spine of an elderly patient with typical deformations.



Geriatric LOR Kit set of 2
LOR inserts
with 2 spinal tubes
XP61-004 – 1020629



Single Geriatric Insert 1
LOR insert
with 1 spinal tube
XP61-003 – 1020800

Delivery contents:

- Base body
- Vertebral bar
- Standard patient LOR insert (loss-of-resistance)
- Skin
- Spinal canal with connector
- Water collection container
- 10 cm tube
- 2 x 100 ml injections with cap
- Extension tube for the hanging-drop technique

Options and replacements:

- Replacement skin (1017892)
- Replacement LOR Kit (1017893)
(2 standard patient LOR inserts with 2 spinal tubes)

For further information, explanatory videos and product manuals visit 3bscientific.com.



3B Scientific GmbH • Rudorffweg 8 • 21031 Hamburg • Germany • 3bscientific.com

Phone: +49 40 73966-0 • Fax: +49 40 73966-100 • E-mail: info@3bscientific.com